


High Pressure Machine Live μ

Technical sheet

- Width : 995,5 mm
- Depth : 1064 mm
- Height : 1855 mm
- Weight : 450 kg
- Voltage : 380Vtri+N+T
- Installed power: 7000W
- Dry air network at 6 bars
- Capacity of the Dewar in LN₂: 9 liters
- Machine connected to a nitrogen container of 60 to 100 litres
- Consumption: 25 liters per day, including cooling
- Alcohol reservoir: 0,5 liter



Révisé		Date	Approuvé par
Echelle		01/06/2017	Jerome
1:5.67	1 / 1		
Montage tolérances Lμ			
TECH		Document n°	